



# REAL-TIME STREAMING ANALYTICS PIPELINE IN AWS FOR HEALTH MONITORING

---

Amit Duwal (080MSDSA004)

Bishwas Dangol (080MSDSA010)

# Introduction

- Real-time streaming analytics allows for continuous processing of data streams, enabling timely decision-making and insights.
- Leverage AWS services to build a scalable, secure, and efficient analytics pipeline for monitoring real-time health data.


# Introduction

- **Why It Matters?:**

- Immediate detection of health anomalies.
- Enhance healthcare outcomes.

- **Key Focus Area:**

- AWS Services (Amazon Kinesis, Amazon DynamoDB, SNS, Lambda, S3)
- Data ingestion, processing, and analysis.
- Real time alert and dashboards(Amazon Open Search)



# Methodology

- **Data Ingestion**
- **Amazon Kinesis** is used to **capture real-time health data** streams. It enables continuous data ingestion from sources like IoT devices (e.g., medical sensors).
- **How it works:**
- Kinesis ingests a high volume of data in real time, splitting the data into shards for efficient processing.
- **Provisioned mode** can be configured with **1 or more shards**, allowing scalable data capture and processing.

# Methodology

- **Cold Storage**
- **Amazon Kinesis Firehose** delivers raw data from **Kinesis Stream** to **Amazon S3** for long-term storage.
- Data from Kinesis Stream is compressed in **GZIP** format using Firehose and stored in S3 for later analysis.
- **Amazon S3** stores large amounts of data with **high durability and availability**, making it ideal for archiving raw health data.
- **Scalable** storage with integrated **lifecycle policies** to manage data cost-effectively.

# Methodology

- **Hot Data**
- **Kinesis Data Analytics** processes real-time data for anomalies (e.g., seizure detection) using SQL-based stream analytics.
- **AWS Lambda** is triggered when a threshold is breached. The Lambda function:
  - Sends real-time alerts to stakeholders via **Amazon SNS**.
  - Writes anomaly details to **Amazon DynamoDB** for further processing and analysis.
- **DynamoDB Table**: Stores patient information, including event timestamps and details of seizures.

# Results And Outcomes

- **Data Storage and Management:**

| REAL TIME SEIZURES TABLE |                 |                               |           |            |                               |                     |
|--------------------------|-----------------|-------------------------------|-----------|------------|-------------------------------|---------------------|
| age                      | currentLocation | eventTimestamp                | heartRate | hospitalId | inspectedAt                   | name                |
| 82                       | null            | 2024-08-19<br>12:50:43.178734 | 161       | null       | 2024-08-19<br>07:05:44.649295 | James<br>Thompson   |
| 61                       | null            | 2024-08-19<br>12:50:27.975594 | 80        | null       | 2024-08-19<br>07:05:29.679300 | Tina<br>Johnson     |
| 67                       | null            | 2024-08-19<br>12:50:26.383464 | 76        | null       | 2024-08-19<br>07:05:28.668891 | Stephen<br>Fuller   |
| 74                       | null            | 2024-08-19<br>12:49:58.911034 | 161       | null       | 2024-08-19<br>07:05:00.665144 | Christian<br>Foster |
| 85                       | null            | 2024-08-19<br>12:49:51.711964 | 150       | null       | 2024-08-19<br>07:04:53.632263 | Bryan<br>Frank      |
| 82                       | null            | 2024-08-19<br>12:50:43.178734 | 161       | null       | 2024-08-19<br>07:05:44.649295 | James<br>Thompson   |



# Methodology

- **Amazon Simple Notification Service (SNS)** sends real-time alerts (e.g., email or SMS) to healthcare providers when critical health events like seizures are detected.
- **How it works:**
  - SNS subscribes healthcare personnel to specific alerts.
  - When a **Lambda function** triggers SNS, it broadcasts notifications to all subscribers.
  - It ensures that **real-time notifications** are delivered instantly, improving response times in emergencies.

# Results And Outcomes

- **Real-Time Notification:**

AWS Notification Message  Inbox x



**Seizure Notification** <no-reply@sns.amazonaws.com>  
to me ▾

 Seizure Alert 

Patient with seizure detected!

Patient Name: Samantha Allen  
Patient ID: eb6d1bf2-ef4e-4a06-82d9-1fe1da9e2305  
Age: 55  
Heart Rate: 107 BPM  
Respiratory Rate: 18 breaths/min  
Oxygen Saturation: 92%  
Seizure Severity: Severe  
Seizure Duration: 178 seconds  
Location: null  
Event Timestamp: 2024-08-19 12:44:35.320600

Please take immediate action.


12:44 PM (5 hours ago)    

# Methodology

- **Amazon OpenSearch Service** is used to **visualize real-time health data** on customizable dashboards.
- **How it works:**
- Data is indexed in **OpenSearch** via **Lambda functions**, enabling dynamic querying and visualization.
- Dashboards are created for monitoring **patient vitals, anomalies**, and trends like seizure occurrence.
- It allows powerful **visualization tools** (e.g., bar charts, pie charts) to display complex datasets, making it easy for healthcare teams to monitor patient health.

# Results And Outcomes

- **Data Visualization Using OpenSearch:**




# Methodology

- **AWS Cloud Development Kit (CDK)** allows defining cloud infrastructure using code.
- It simplifies **cloud deployment** and offers a **high-level abstraction** using familiar programming languages.
- **Benefits:**
  - **Infrastructure as Code (IaC)** for easier **replication** and scalability.
  - Integration with multiple AWS services.
  - **Automated** health monitoring system setup.
  - **Customizable pipelines** for different health monitoring use cases.


# Methodology

```
102 // =====
103 // OpenSearch Module:
104
105 const opensearchDomain = new Domain(this, 'Domain', {
106   version: EngineVersion.OPENSEARCH_2_13,
107   domainName: 'opensearch-domain',
108   capacity: {
109     dataNodes: 1,
110     dataNodeInstanceType: 't3.medium.search',
111     multiAzWithStandbyEnabled: false,
112   },
113   ebs: {
114     volumeSize: 20,
115   },
116   enforceHttps: true,
117   nodeToNodeEncryption: true,
118   encryptionAtRest: {
119     enabled: true,
120   },
121   fineGrainedAccessControl: {
122     masterUserName: process.env.OPENSEARCH_USER,
123     masterUserPassword: cdk.SecretValue.unsafePlainText(process.env.OPENSEARCH_PASSWORD as string),
124   },
125   accessPolicies: [
126     new iam.PolicyStatement({
127       actions: ['es:*'],
128       resources: ['*'],
129       effect: iam.Effect.ALLOW,
130       principals: [new iam.AnyPrincipal()],
131     )],
132   removalPolicy: cdk.RemovalPolicy.DESTROY,
133 });
134
```

# Results And Outcomes




# Results And Outcomes







# Results And Outcomes



# Results And Outcomes




# Results And Outcomes




# Results And Outcomes

- **Dashboard Creation:**



# Results And Outcomes

- Dashboard Creation:



# DEMO

# Conclusion

- The project provides **Continuous real-time data ingestion** from wearable devices and medical sensors along with **Timely alerts** and notifications via SNS for rapid response to health changes. Also provides **Secure and scalable storage** in S3 for both real-time and historical data and OpenSearch for **Visualization**.
- Enhances **patient outcomes** and supports **personalized treatment plans**. **Optimizes healthcare resources** and improves overall healthcare efficiency.

Thank You